

**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. CNP- WB-548**

# Anti PDCD4

**Product type** Primary antibodies

**Immunogen** PDCD4 (AAH26104.1, 1 a.a. ~ 469 a.a) full-length human protein.

**Form** Liquid  
Protein A purified

**Volume** 100ul

**Specificity** Human


**Host** Rabbit

**Immunogen Sequence** MDVENEQILNVNPADPDNLSDSLFSGDEENAGTEEVKNEINGNWISAYSINEARINAKAK  
RRLRNSSRDSDGRGDSVSDSGSDALRSGLTVPTSPKGRLLDRRSRSGKGRGLPKKGGAGG  
KGVWGTPGQVYDVEEVDVKDPNYDDDQENCYETVVLPLDERAFEKTLTPIIQEYFEHGD  
TNEVAEMLRDLNLGEMKSGVPVLAVSLALEGKASHREMTSKLLSDLCGTVMSTTDVEKSF  
DKLLKDLPELALDTPRAPQLVGQFIARAVGDGILCNTYIDSYKGTVDVCVQARAALDKATV  
LLSMSKGGKRKDSVWGSGGGQQSVNHLVKEIDMLLKEYLLSGDISEAEHCLKELEVPHFH  
HELVYEAIIIMVLESTGESTFKMILDLLKSLWKSSTITVDQMKGRIYNEIPDINLDVP  
HSYSVLERFVEECFQAGIISKQLRDLCP SRGRKRFVSEGDGGRLKPESY

**Gene Alias** H731, MGC33046, MGC33047

**Application notes** **WB**  
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~51.80 kDa using PDCD4 transfected 293T lysate.

**Western Blotting (Transfected lysate)**



Western Blot analysis of PDCD4 expression in transfected 293T cell line by PDCD4 rabbit polyclonal antibody.

Lane 1: PDCD4 transfected lysate(51.80 KDa).

Lane 2: Non-transfected lysate.

**Storage**

For continuous use, store at 2-8°C or one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619