

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-595

Anti IFNG

Product type Primary antibodies

Immunogen IFNG (AAH70256.1, 1 a.a. ~ 166 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MKYTSYILAFQLCIVLGS LGCYCQDPYVKEAENLKKYFNAGHSDVADNGTLFLGILKNWK
EESDRKIMQSQIVSFYFKLFKNFKDDQSIQKSVETIKEDMNVKFFNSNKKKRDDFEKLTN
YSVTDLNVQRKAIHELIQVMAELSPA AKTGKRKRSQMLFRGRRASQ


Gene Alias IFG, IFI

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~19.30 kDa using IFNG transfected 293T lysate..

Western Blotting (Transfected lysate)



Western Blot analysis of IFNG expression in transfected 293T cell line by IFNG rabbit polyclonal antibody.

Lane 1: IFNG transfected lysate(19.30 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619