

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-596

Anti HFE2

Product type Primary antibodies

Immunogen HFE2 (NP_660320.3, 1 a.a. ~ 313 a.a) full-length human protein

Form Liquid

Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MIQHNC SRQGPTAPPPPRGPALPGAGSGLPAPDPCDYEGRFSRLHGRPPGFLHCASFGDP
HVRSFHHHFHTCRVQGA WPLLDNDFLFVQATSSPMALGANATATRKLTIIIFKNMQECIDQ
KYYQAEVDNLPVAFEDGSINGGDRPGGSSLSIQTANPGNHVEIQAAAYIGTTIIIRQTAGQ
LSFSIKVAEDVAMAFSAEQDLQLCVGGCPPSQRLSRSENRNRGAITIDTARRLCKEGLPV
EDAYFHSCVFDVLISGDPNFTVAAQALEDARAFLPDLEKLHLFSPDAGVPLSSATLLAP
LLSGLFVLWLCIQ

Gene Alias HFE2A, HJV, JH, MGC23953, RGMC


Application

notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~33.70 kDa using HFE2 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of HFE2 expression in transfected 293T cell line by HFE2 rabbit polyclonal antibody.

Lane 1: HFE2 transfected lysate(33.70 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of HFE2 expression in mouse spleen (Predicted Mw: 44.8480 kDa).

Western Blotting (Tissue lysate)



HFE2 rabbit polyclonal antibody. Western Blot analysis of HFE2 expression in mouse spleen.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619