



POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-WB-600

Anti NUDT5

Product type	Primary antibodies
Immunogen	NUDT5 (NP_054861.2, 1 a.a. ~ 219 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MESQEPTESSQNGKQYIISEELISEGKWVKLEKTTYMDPTGKTRTWESVKRTRRKEQTAD GVAVIPVLQRTLHYECIVLVKQFRPPMGGYCIEFPAGLIDDGETPEAAALRELEEETGYK GDIAECSPAVCM DPGLSNCTIHIVTVTINGDDAENARPKPKPGDGEFVEVISLPKNDLLQ RLDALVAEEHLTVDARVYSYALALKHANAKPFEVPFLKF
Gene Alias	YSA1, YSA1H, hYSAH1
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~24.30 kDa using NUDT5 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of NUDT5 expression in transfected 293T cell line by NUDT5 rabbit polyclonal antibody.

Lane 1: NUDT5 transfected lysate(24.30 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619