

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-WTIP-009

Anti WTIP

Product type	Primary antibodies
Immunogen	WTIP (AAI60027.1, 1 a.a. ~ 430 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

**Immunogen
Sequence**


MQRSRAGADEAALLLAGLALRELEPGCGSPGRGRRGPRPGPGDEAAPALGRRGKGS GGPE
AGADGLSRGERGPRRAAVPELSAQPAGSPRASLAGSDGGGGGGGSARSSGISLG YDQRHGS
PRSGRSDPRPGPGPPSVGSARSSVSSLGSRGSAGAYADFLPPGACPAPARSPEPAGPAPF
PLPALPLPPGREGGPSAAERRLEALTRELERALEARTARDYFGICIKCGLGIYGAQQACQ
AMGSLYHTDCFTCDSCGRRLRGKAFYNVGEKVYCQEDFLYSGFQQTADKCSVC GHLIMEM
ILQALGKSYHPGCFRCSVCNECLDGVPFTVDVENNIYCVRDYHTVFAPKCASCARPILPA
QGCETTIRVVSMRDRDYHVACYHCEDCGLQLSGEEGRRCYPLAGHLLCRRCHL RRLQPGPL
PSPTVHVTEL

**Application
notes**

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~47.30 kDa using WTIP transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of WTIP expression in transfected 293T cell line by WTIP rabbit polyclonal antibody.

Lane 1: WTIP transfected lysate(47.30 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619