

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-ZFYVE9-146

Anti ZFYVE9

Product type	Primary antibodies
Immunogen	ZFYVE9 (NP_015562.1, 1 a.a. ~ 762 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit


Immunogen Sequence

MENYFQAEAYNLDKVLDEFEQNEDETVSSTLLDTKWNKILDPPSHRLSFNPTLASVNESA
 VSNESQPQLKVFSLAHSAPLTTEEDHCANGQDCNLNPEIATMWIDENAVAEDQLIKRNY
 SWDDQCSAVEVGEKKCGNLACLPDEKNVLVAVMHNC DKRTLQNDLQDCNNYNSQSLMDA
 FSCSLDNENRQTDQFSFSINESTEKDMNSEKQMDPLNRPKTEGRSVNHLCPTSSDSLAV
 CSPSQLKDDGSIGRDPMSAITS LTVDSVISSQGTGDCPAVKKQENYIPDEDLTGKISSP
 RTDLGSPNSFSHMSEGILMKKEPAEESTTEESLRSGPLLLLKPDMPNGSGRNND CERCSD
 CLVPNEVRADENEGYEHEETLGTTEFLNMTEHFSESQDMTNWKLTKLNEMNDSQVNEEKE
 KFLQISQPEDTNGDSGGQCVGLADAGLDLKGTCISESEECDFSTVIDTPAANYLSNGCDS
 YGMQDPGVSVFPKTLPSKEDSVTEEKEIEESKSECYSNIYEQRGNEATEGSGLLLNSTGD
 LMKKNYLHNFC SQVPSVLGQSSPKVVASLPSISVPFGGARPKQPSNLKLQIPKPLSDHLQ
 NDFPANSNGNNTKNKNDILGKAKLGENSATNVCSPSLGNISNVD TNGEHLESYEA EISTRP
 CLALAPDSPDNDLRAGQFGISARKPFTTLGEVAPVWVPDSQAPNCMKCEARFTFTKRRHH
 CRACGKVFCASCCSLKCKLLYMDRKEARVCVICHSVLMNGKY

Gene Alias MADHIP, NSP, SARA, SMADIP

Application notes **WB**
 Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~83.70 kDa using ZFYVE9 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of ZFYVE9 expression in transfected 293T cell line by ZFYVE9 rabbit polyclonal antibody.

Lane 1: ZFYVE9 transfected lysate(83.70 KDa).
 Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of ZFYVE9 transfected lysate.

Immunoprecipitation



Immunoprecipitation of ZFYVE9 transfected lysate using ZFYVE9 rabbit polyclonal antibody and Protein A Magnetic Bead, and Immunoblotted with ZFYVE9 rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619