



POLYCLONAL ANTIBODY

For research use only, Not for diagnostic use.

Catalog No. EXO-AB-01

Anti Bovine Milk Exosome

Product type	Primary antibody
Immunogen	Bovine
Rased in	Rabbit
Myeloma	-
Clone number	-
Isotype	-
Source	Serum
Purification	Ammonium sulfate precipitation Milk Whey protein absorbed
Buffer	PBS(-) 15ppm Proclin300
Concentration	1 mg / mL
Volume	100 ul
Label	Unlabeled
Specificity	Whole Exosome of bovine Milk
Cross reactivity	Bovine
Storage	Store below - 20°C Aliquot to avoid cycles of freeze/thaw.

Application notes • **Western blotting:** 1/1000

Recommended dilutions

References -

RELATED PRODUCT:

Product Name	Quantity	Maker	Cat#
Milk exosome ,Bovine	1 SET (100uLx10 vial)	CSR	EXBM100L
Milk exosome ,Bovine	1 SET (1mLx10 vial)	CSR	EXBM1000L
Milk Exosome ELISA Kit ,Bovine	1 KIT (96test)	CSR	EXBMEL

For research use only, Not for diagnostic use.



COSMO BIO Co., LTD.

[JAPAN]
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

[Outside Japan]
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com