

## MONOCLONAL ANTIBODY

*For research use only, Not for diagnostic use.***Catalog No. SHI-EXO-M02-TF2****Anti CD63 [ Clone : 8A12]****BACKGROUND**

CD63 (also known as LAMP-3, Melanoma-associated antigen ME491, TSPAN30, MLA1 and OMA81H) is a cell surface glycoprotein which belongs to the tetraspanin superfamily. CD63 is known to complex with integrins. CD63 is expressed on activated platelets, monocytes and macrophages, and is weakly expressed on granulocytes, T cell and B cells.

<b>Product type</b>	Primary antibody
<b>Immunogen</b>	Human CD63 (a.a. 104-202)
<b>Raised in</b>	Mouse
<b>Myeloma</b>	-
<b>Clone number</b>	8A12
<b>Isotype</b>	IgG2a
<b>Source</b>	Culture supernatant
<b>Purification</b>	Affinity purified by Protein G
<b>Buffer</b>	Phosphate buffered saline containing 15 ppm ProClin
<b>Concentration</b>	1 mg / mL
<b>Volume</b>	100 ul
<b>Label</b>	Tide Fluor™ 2WS Ex:502nm, Em :525nm
<b>Specificity</b>	Human CD63
<b>Reactivity</b>	Human. No cross-reactivity with Rat and Mouse; Other species are not tested.
<b>Storage</b>	Store below - 40°C Aliquot to avoid cycles of freeze/thaw.

**RELATED PRODUCT**

Product Name	Application	Quantity	Maker	Cat#
Anti CD9 (12A12) Monoclonal Antibody	WB, IP, ELISA, FC, EM, IF	100 ug	CAC	SHI-EXO-M01
Anti CD63 (8A12) Monoclonal Antibody	WB, IP, ELISA, FC, EM, IF	100 ug	CAC	SHI-EXO-M02
Anti CD81 (12C4) Monoclonal Antibody	WB, IP, ELISA, FC, EM, IF	100 ug	CAC	SHI-EXO-M03

*For research use only. Not for clinical diagnosis.***COSMO BIO Co., LTD.**

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>

**COSMO BIO USA**

【Outside Japan】

2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobiousa.com](mailto:info@cosmobiousa.com)  
Phone/FAX: (+1) 760-431-4600  
URL: [www.cosmobiousa.com](http://www.cosmobiousa.com)