



MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. SU-IZ-M03

Anti human IL-13Ra D1 domain (SS4B)

BACKGROUND

Anti-SCCA2 antibody is a rat monoclonal antibody which obtained from the immunization with purified, E. coli-derived, recombinant human SCCA2. This antibody can be used for the detection of human SCCA2 by immunoprecipitation and ELISA with no cross-reaction to SCCA1.

Product type	Primary antibodies
Host	Rat
Source	Mouse ascites
Form	Liquid
	Purified IgG with 0.05% NaN ₃ as a preservative.
Volume	100 ul
Concentration	0.47 mg/ml
Specificity	IL-13 receptor α 1D1 domain
Antigen	Recombinant human D1 domain expressed by E. coli
Clone	SS4B
Isotype	IgG2a/k

Application notes	FC
	Recommended dilutions
	FC, 1-10 ug/ml
	Optimal dilutions/concentrations should be determined by the end user.
Cross reactivity	Human
Storage	Store below -20°C (below -70°C for prolonged storage). Aliquot to avoid cycles of freeze/thaw.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618 FAX : +81-3-5632-9619