



MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. SU-IZ-M07

Anti human SCCA1/SCCA2

BACKGROUND

Anti-SCCA1/2 antibody is a rat monoclonal antibody which obtained from the immunization with purified, E. coli-derived, recombinant human SCCA1. This antibody can be used for the detection of human SCCA1 and SCCA2 by immunoprecipitation and ELISA.

Product type	Primary antibodies
Host	Rat
Source	Culture supernatant
Form	Liquid Purified IgG with 0.05% NaN ₃ as a preservative.
Volume	100 ul
Concentration	0.3 mg/l
Specificity	SCCA1, SCCA2
Antigen	Recombinant SCCA1 expressed by E. coli.
Clone	SS6A
Isotype	IgG2b/k

Application notes	IP, ELISA Recommended dilutions IP, 1-10 ug/ml ELISA, 1-10 ug/ml Optimal dilutions/concentrations should be determined by the end user.
Cross reactivity	Human
Storage	Store below -20°C (below -70°C for prolonged storage). Aliquot to avoid cycles of freeze/thaw.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618 FAX : +81-3-5632-9619