



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. SU-IZ-P01

Anti EGFP

BACKGROUND

Anti EGFP antibody is a rabbit polyclonal antibody which obtained from the immunization with purified, E. coli-derived, recombinant EGFP. This antibody has been tested and confirmed for use in immunoblotting against enhanced GFP.

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid Serum with 0.05% NaN ₃ as a preservative.
Volume	100 ul
Concentration	
Specificity	GFP
Antigen	Recombinant EGFP expressed by E. coli.

Application notes	WB Recommended dilutions Western blotting, 1/1,000 to 1/5,000 Optimal dilutions/concentrations should be determined by the end user.
Storage	Store below -20°C (below -70°C for prolonged storage e). Aliquot to avoid cycles of freeze/thaw.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618 FAX : +81-3-5632-9619