



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. SU-IZ-P04

Anti human SCCA

BACKGROUND

Anti-SCCA antibody is a rabbit polyclonal antibody which obtained from the immunization with purified, E. coli-derived, recombinant human SCCA2. This antibody can be used for the detection of SCCA by immunoblotting and immunostaining.

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid Serum with 0.05% NaN ₃ as a preservative.
Volume	100 ul
Specificity	SCCA1, SCCA2
Antigen	Recombinant SCCA2 expressed by E. coli.

Application notes	WB, IHC Recommended dilutions Western blotting, 1/100 to 1/1,000 Immunohistochemistry, 1/50 to 1/500 Optimal dilutions/concentrations should be determined by the end user.
Cross reactivity	Human
Storage	Store below -20°C (below -70°C for prolonged storage). Aliquot to avoid cycles of freeze/thaw.
References	1) Mitsuishi, K. et al., Clin. Exp. Allergy. 35:1327-1333, 2005