



MONOCLONAL ANTIBODY

Catalog No. TNL-001-CAM

Anti *Drosophila* CaMKII (dCaMKII)

Product type	Primary antibodies
Host	Mouse
Source	Ascites
Form	Liquid containing 0.02% NaN ₃ as preservative
Volume	100 µl
Concentration	
Specificity	
Antigen	Recombinant purified His-tagged <i>Drosophila</i> CaMKII (490 amino acid isoform)
Clone	18
Isotype	IgG2b

Application notes WB, IHC

Recommended use

Recommended dilutions

Western blotting, 1/1,000 to 1/4,000 Predicted molecular weight: 55-, 58-, 60-kDa
Immunohistochemistry, 1/1,000 to 1/4,000

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage

Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

References

- 1) Takamatsu, Y., Kishimoto, Y., and Ohsako, S. Immunohistochemical study of Ca²⁺/calmodulin-dependent protein kinase II in the *Drosophila* brain using a specific monoclonal antibody. Brain Res. 974:99-116. 2003

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

【JAPAN】
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com