

**Catalog No.CBX00395**

## Mouse monoclonal antibody

## Anti-Human LSM8

## ■ Formulation

Mouse monoclonal anti-human **LSM8** antibody in PBS (3.0 mM KCl, 1.5 mM  $\text{KH}_2\text{PO}_4$ , 140 mM NaCl, 8.0 mM  $\text{Na}_2\text{HPO}_4$  (pH 7.4)) containing 1% bovine serum albumin (BSA) and 0.05% sodium azide ( $\text{NaN}_3$ ).

### ■ Antibody concentration

100µg/ml

## ■ Storage

Store at 2-8°C for up to one year.  
We recommend storing at -20°C for long-term storage.  
Avoid repeat freezing and thawing cycles.

## ■ Preparation

This antibody was purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum.

## ■ Sterility

0.22μm membrane

## ■ Applications

WB, Dot Blot

## ■ Disposal

This antibody solution contains sodium azide ( $\text{NaN}_3$ ) as a preservative. There is a potential hazard that  $\text{NaN}_3$  reacts with copper or lead to produce an explosive compound. For safe disposal, the vial has to be washed thoroughly with water.

## ■ Safety warnings and precautions

Caution must be taken to avoid contact with skin or eyes. In such a case, rinse thoroughly at once with water. Do not ingest, inhale, or swallow. Seek medical attention immediately.

Wear appropriate protective clothing such as laboratory overalls, safety glasses and gloves.

It is strongly advised that this product should be handled by people who have been well trained in

laboratory techniques and that it is handled with care pursuant to the principles of good laboratory practice. All chemicals are deemed potentially harmful.

The vial is prone to fall over. Use caution, especially when the lid is off.

**Lot No. 2527C6a-1**


**Clone No.** 2527C6a

**Antibody class : IgG1**

## Immunogen : Recombinant



Western blot analysis of immunized recombinant protein, using anti- LSM8 monoclonal antibody.

N  C

Arrow indicates the region of immunized recombinant protein carrying 50-200 amino acids.

FOR RESEARCH USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES.

**manufactured by BMR**