

Clonogene LLC

Product Name Dicer antibody (13D6) / **CG-006**

Product type Primary antibody

Description Mouse monoclonal to Dicer

Target protein Dicer

Immunogen Synthetic peptide from N-terminus of Dicer

Specificity reacts with at least human, mouse and rat Dicer.

Raised in mouse

Clonality monoclonal

Isotope IgG1

Purity Hybridoma supernatant

Storage buffer N/A

Form Liquid

Tested applications: Western blot. Other apps pending

Application notes (see key)

WB: Use at a dilution of 0.5ug/mL (1/100 dilution of sup.) in 1x BLOTTO or 3% BSA in PBS (with ECL or ECL plus) Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Positive control Detects 10ng of recombinant Dicer (using ECL or ECL Plus) in under one minute of exposure to film.

Relevance Excellent in western blots and potentially useful in other applications

References

N/A just released