



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. COP-080051

Anti Yeast Sis1 (DnaJ)

Product type	Primary antibodies
Host	Rabbit
Source	Serum
Form	Liquid
	0.1% NaN ₃ as a preservative.
Volume	100 µl
Concentration	
Specificity	
Antigen	synthetic peptide: NH ₂ - C+SLNDAQKRAIDENF -COOH
Isotype	

Application notes WB, IP

Recommended use

Recommended dilutions

Western blotting: ~1/5,000 Predicted molecular weight: 38kDa


Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage Store below -30°C (below -80°C for prolonged storage).
Aliquot to avoid cycles of freeze/thaw.

References



For research use only. Not for clinical diagnosis.

Manufactured by Operon Biotechnologies, K.K



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619