



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CZA-002

Anti KIAA1199


Product type	Primary antibodies
Immunogen	KIAA1199 (AAH20256.1, 1 a.a. ~ 992 a.a) full-length human protein.
Form	Liquid Protein G purified
Quantity	100ug
Buffer	PBS
Preservative	ProClin (15ppm)
Specificity	Human
Host	Rabbit

Immunogen Sequence	MGAAGRQDFLFKAMLTISWLTLCFPGATSTVAAGCPDQSPQLQPWNPGHDQDHHVHIGQ GKTLTLLTSSATVYSIHISEGGKLVKDHDEPIVLRTRHILIDNGGELHAGSALCPFQGNF TIILYGRADEGIQPDYGLKYIGVGKGGALELHGQKLSWTFLNKTLHPGGMAEGGYFF ERSWGHARGVIVHVIDPKSGTVIHSDFDFTYRSKKESERLVQYLNAVDPGRILSVAVNDEG SRNLDDMARKAMTKLGSKHFLHLGFRHPWSFLT VKGNPSSSVEDHIEYHGHRGSAAARVF KLFQTEHGGEYFNVSLSSEWVQDVEWTEWFDHDKVSQTKGGEKISDLWKAHPGKICNRPID IQATTMDGVNLSDEVVYKKGQDYRFACYDRGRACRSYRVRFLCGKPVPRPKLTVTIDTNVN STILNLEDNVQSWKPGD TLVIASTDYSMYQAEFQVLP CRSCAPNQVKVAGKPMYHIGE EIDGVDMRAEVGLLSRNIIVMGEMEDKCYPYRNHICNFFDFDTFGGHIKFALGFKAAHLE GTELKHMGGQLVGGYPIHFHLAGDVDERGGYDPPTYIRDLSIHHTFSRCVTVHGSNGLLI KDVVGYNSLGHCFFTEDGPEERNTFDHCLG LLVKSGTLLPSDRDSKMCKMITEDSYPGYI PKPRQDCNAVSTFWMANPNNNLINCAAAGSEETGFWFIFHHVPTGPSVGMYSYSPGYSEHIP LGKFYNNRAHSNYRAGMIIDNGVKTTTEASAKDKRPFLSIISARYSPHQDADPLKPREPAI IRHFIAYKNQDHGAWLRGGDVWLDSCR FADNGIGLTLASGGTFPYDDGSKQEIKNSLFVG ESGNVGTMMDNRIWGPGLDHSGR TLP I GQNFPIRGIQLYDGPINIQNCTFRKFVALEG RHTSALAFRLNNAWQSCPHNNVTGIAFEDVPITSRVFFGEPGPWFNQLDMDGDKTSVFHD VDGSVSEYPGSYLTKNDKWHSLASKAASGPSG
---------------------------	--

Gene Alias TMEM2L

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~110.40 kDa using KIAA1199 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of KIAA1199 expression in transfected 293T cell line by KIAA1199 rabbit polyclonal antibody.


Lane 1: KIAA1199 transfected lysate(110.40 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of KIAA1199 transfected lysate.

Immunoprecipitation



Immunoprecipitation of KIAA1199 transfected lysate using KIAA1199 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with KIAA1199 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

Production

Cosmo Bio Co., Ltd., Japan

For research use only, Not for diagnostic use.



COSMO BIO Co., LTD.

[JAPAN]
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU, TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

[Outside Japan]
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com