



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CZA-003

Anti TP53

Product type	Primary antibodies
Immunogen	TP53 (NP_000537.2, 1 a.a. ~ 393 a.a) full-length human protein.
Form	Liquid Protein G purified
Quantity	100ug
Buffer	PBS
Preservative	ProClin (15ppm)
Specificity	Human
Host	Rabbit

Immunogen Sequence
MEEPQSDPSVEPPLSQETFSDLWKLLPENNVLSPLPSQAMDDLMLSPDDIEQWFTEDPGP
DEAPRMPEAAPRVAPAPAAPTPAAPAPAPSWPLSSSVPSQKTYQGSYGFRLGFLHSGTAK
SVTCTYSPALNKMFCQLAKTCPVQLWVDSTPPPGRVTRAMAIYKQSQHMTEVVRRCPHHE
RCSDSDGLAPPQHLIRVEGNLRVEYLDDRNTFRHSVVVPYEPPEVGS DCTTIHYNMCNS
SCMGGMNRRPILTIITLEDSSGNLLGRNSFEVRVCACPGRRRTEENLRKKGEPHHELP
PGSTKRALPNNTSSSPQPKKKPLDGEYFTLQIRGRERFEMFRELNEALELKDAQAGKEPG
GSRAHSSHLSKSKKGQSTSRHKKLMFKTEGPDS


Gene Alias FLJ92943, LFS1, TRP53, p53

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~43.23 kDa using TP53 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of TP53 expression in transfected 293T cell line by TP53 rabbit polyclonal antibody.

Lane 1: TP53 transfected lysate(43.23 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

Production

Cosmo Bio Co., Ltd., Japan

For research use only, Not for diagnostic use.



COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU, TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com