



**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. CZA-005**

**Anti PH4**

<b>Product type</b>	Primary antibodies
<b>Immunogen</b>	P4HTM (NP_808807.2, 1 a.a. ~ 563 a.a) full-length human protein.
<b>Form</b>	Liquid Protein G purified
<b>Quantity</b>	100ug
<b>Buffer</b>	PBS
<b>Preservative</b>	ProClin (15ppm)
<b>Specificity</b>	Human
<b>Host</b>	Rabbit

**Immunogen Sequence**


MAAAVTGQRPETAAAEASRPQWAPPDHCQAQAAAGLGDGEDAPVRPLCKPRGICSRAY  
FLVLMVFVHLYLGNVLALLLFVHYSNGDESSDPGPQHRAQGGPEPTLGPLTRLEGIKVG  
HERKVQLVTRDRDHFIRTLCLKPLLEIPGFLTDEECRLIIHLAQMKGLQRSQILPTEEYE  
EAMSTMQVSQDLFRLLDQNRDGHLLQREVLAQTRLGNGWWMTPESIQEMYAAIKADPDG  
DGVLSLQEFNSMDLRDFHKYMRSHKAESSELVRNSHHTWLYQGEGAHHIMRAIRQRVLR  
TRLSPEIVELSEPLQVVRYGEGGHYHAHVDSGPVYPETICSHTKLVANESVPFETSCRQV  
SPNWGLPSILRPGTPMTQAQPCTVGVPLGMGPGDHWVIVSPWEHPQLGTCSVPPLPYSY  
MTVLFYLNVTGGGETVFPVADNRTYDEMSLIQDDVDLRDTRRHCDKGNLRVKPQQGTAV  
FWYNYLPDGQGWVGDVDDYSLHGGCLVTRGTKWIANNWINVDPSRARQALFQQEMARLAR  
EGGTDSSQPEWALDRAYRDARVEL

**Gene Alias** EGLN4, HIFPH4, P4H-TM, PH-4, PH4, PHD4

**Application notes**

**WB**  
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~63.10 kDa using P4HTM transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of P4HTM expression in transfected 293T cell line by P4HTM rabbit polyclonal antibody.

Lane 1: P4HTM transfected lysate(63.10 KDa).  
Lane 2: Non-transfected lysate.

**Storage** For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours

**Production** Cosmo Bio Co., Ltd., Japan

For research use only, Not for diagnostic use.



**COSMO BIO Co., LTD.**

**【JAPAN】**  
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU, TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>



**COSMO BIO USA**

**【Outside Japan】**  
2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobioussa.com](mailto:info@cosmobioussa.com)  
Phone/FAX: (+1) 760-431-4600  
URL: [www.cosmobioussa.com](http://www.cosmobioussa.com)