



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CZA-008

Anti GALNS


Product type	Primary antibodies
Immunogen	GALNS (NP_000503.1, 1 a.a. ~ 522 a.a) full-length human protein.
Form	Liquid Protein G purified
Quantity	100ug
Buffer	PBS
Preservative	ProClin (15ppm)
Specificity	Human
Host	Rabbit

Immunogen Sequence
 MAAVVAATRWWQLLLVLSAAGMGASGAPQPPNILLLLMDDMGWGD LGVYGEPSRETPNLD
 RMAAEGLLFPNFYSANPLCSPSRAALLTGRLPIRNGFYTTNAHARNAYTPQEIVGGIPDS
 EQLLPELLKKAGYVSKIVGKWHLGHRPQFHPLKHGFDEWFGSPNCHFGPYDNKARNIPV
 YRDWEMVGRYYEEFPINLKTGEANLTQIYLQEALDFIKRQARHHPFFLYWAVDATHAPVY
 ASKPFLGTSQRGRYGDVREIDDSIGKILELLQDLHVADNTFVFFTS DNGAALISAPEQG
 GSNPFLCGKQTTFEGGMREPALAWWPGHV TAGQVSHQLGSIMDLFTTSLALAGLTPPSD
 RAIDGLNLLPTLLQGRLMDRPIFYR GDTLMAATLGQHKAHFWTWTNSWENFRQGIDFCP
 GQNVSGVTTHNLEDHTKLPLIFHLGRDPGERFPLSFASAEYQEALSRITSVVQQHQEALV
 PAQPQLNVCNWAVMNWAPPGCEKLGKCLTPPESIPKKCLWSH

Gene Alias FLJ17434, FLJ42844, FLJ98217, GALNAC6S, GAS, MPS4A

Application notes **WB**
 Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~58.00 kDa using GALNS transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of GALNS expression in transfected 293T cell line by GALNS rabbit polyclonal antibody.

Lane 1: GALNS transfected lysate(58.00 KDa).

Lane 2: Non-transfected lysate.

IP
Immunoprecipitation



Immunoprecipitation of GALNS transfected lysate using GALNS rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with GALNS rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

Production

Cosmo Bio Co., Ltd., Japan

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619