



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CZA-009

Anti LGALS3BP


Product type	Primary antibodies
Immunogen	LGALS3BP (NP_005558.1, 1 a.a. ~ 585 a.a) full-length human protein
Form	Liquid Protein A purified
Quantity	100ug
Buffer	PBS
Preservative	ProClin (15ppm)
Specificity	Human
Host	Rabbit

Immunogen Sequence
 MTPPRLFWVWLLVAGTQGVNDGDMRLADGGATNQGRVEIFYRGQWGTVCNDLWDLTDASV
 VCRALGFENATQALGRAAFGQGGSPIMLDEVQCTGTEASLADCKSLGWLKSNCRHERDAG
 VVCTNETRSTHTLDLSRESEALGQIFDSQRGCDLSISVNVQGEDALGFCGHTVILTANL
 EAQALWKEPGSNVTMSVDAECVPMVRDLLRYFYSTRIDITLSSVKCFHKLASAYGARQLQ
 GYCASLFAILLPQDPSFQMPLDLYAYAVATGDALLEKLCLQFLAWNFEALTQAEAWPSVP
 TDLLQLLLPRSDLAVPSELALLKAVDTWSWGERASHEEVEGLVEKIRFPMMLPEELFELQ
 FNLSLYWSHEALFQKKTLLQALEFHTVPFQLLARYKGLNLTEDTYKPRIYTSPTWSAFVTD
 SWSARKSQLVYQSRGRLVKYSSDYFQAPSDYRYYPYQSFQTPQHPSFLFQDKRVSWSL
 VYLPTIQSCWNYGFSCSSDELPLVGLTKSGGSDRTIAYENKALMLCEGLFVADVTDFEGW
 KAAIPSALDTNSSKSTSSFPCPAGHFNGFRTVIRPFYLTNSSGVD

Gene Alias 90K, BTBD17B, MAC-2-BP

Application notes **IP**
 Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of LGALS3BP transfected lysate

Immunoprecipitation



Immunoprecipitation of LGALS3BP transfected lysate using LGALS3BP rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with LGALS3BP purified mouse polyclonal antibody.

Storage For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

Production Cosmo Bio Co., Ltd., Japan

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619