



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CZA-009

Anti LGALS3BP


Product type	Primary antibodies
Immunogen	LGALS3BP (NP_005558.1, 1 a.a. ~ 585 a.a) full-length human protein
Form	Liquid Protein A purified
Quantity	100ug
Buffer	PBS
Preservative	ProClin (15ppm)
Specificity	Human
Host	Rabbit

Immunogen Sequence
 MTPPRLFWWLLVAGTQGVNDGDMRLADGGATNQGRVEIFYRGQWGTVCNDLWDLTDASV
 VCRALGFENATQALGRAAFGQGGSPIMLDEVQCTGTEASLADCKSLGWLKSNCRHERDAG
 VVCTNETRSTHTLDLSRELSEALGQIFDSQRGCDLSISVNVQGEDALGFCGHTVILTANL
 EAQALWKEPGSNVTMSVDAECVPMVRDLLRYFYSTRIDITLSSVKCFHKLASAYGARQLQ
 GYCASLFAILLPQDPSFQMPLDLYAYAVATGDALLEKLCLQFLAWNFEALTQAEAWPSVP
 TDLLQLLLPRSDLAVPSELALLKAVDTWSWGERASHEEVEGLVEKIRFPMMLPEELFELQ
 FNLSLYWSHEALFQKKTLLQALEFHTVPFQLLARYKGLNLTEDTYKPRIYTSPTWSAFVTD
 SWSARKSQLVYQSRRGPLVKYSSDYFQAPSDYRYYPYQSFQTPQHPSFLFQDKRVSWSL
 VYLPTIQSCWNYGFSCSSDELPLVGLTKSGGSDRTIAYENKALMLCEGLFVADVTFEGW
 KAAIPSALDTNSSKSTSSFPCPAGHFNGFRTVIRPFYLTNSSGVD

Gene Alias 90K, BTBD17B, MAC-2-BP

Application notes **IP**
 Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of LGALS3BP transfected lysate

Immunoprecipitation



Immunoprecipitation of LGALS3BP transfected lysate using LGALS3BP rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with LGALS3BP purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

Production

Cosmo Bio Co., Ltd., Japan

For research use only, Not for diagnostic use.



COSMO BIO Co., LTD.

【JAPAN】
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU, TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com