

*For test research. It cannot be used for any other purpose.*

**Catalog No. RMP-0002**

## HA Peptide

Peptide name	HA Peptide
Purity	≥95%
Sequence	[H]YPYDVDPYA[OH]
Molecular formula	C <sub>53</sub> H <sub>67</sub> N <sub>9</sub> O <sub>17</sub>
molecular weight	1,102.15 Da.
fluorescence	–
Lot#	PP00000045-002
Counter ion	TFA
Dissolved solvent	Water
Quantity	2.0 mg (1.0 mg×2)

Storage Store below -20 °C. For longterm storage, recommend -80°C or below.  
Please use the dissolved sample in a short time.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

【JAPAN】  
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>



**COSMO BIO USA**

【Outside Japan】  
2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobioussa.com](mailto:info@cosmobioussa.com)  
URL: [www.cosmobioussa.com](http://www.cosmobioussa.com)  
Phone/FAX: (+1) 760-431-4600

試験研究用となります。それ以外の用途のご使用はできません。

Catalog No. RMP-0002

## HA Peptide

ペプチド名	HAペプチド
純度	≥95%
ペプチド配列	[H]YPYDVPDYA[OH]
分子式	C <sub>53</sub> H <sub>67</sub> N <sub>9</sub> O <sub>17</sub>
分子量	1,102.15 Da.
蛍光・標識等	—
Lot#	PP00000045-002
カウンターイオン	トリフルオロ酢酸
推奨溶解溶媒	水
内容量	2.0 mg (1.0 mg×2)

使用・保管方法 -20℃以下で保管してください。長期保管する場合は-80℃以下を推奨します。溶解したサンプルは短時間でご使用下さい。



コスモ・バイオ株式会社  
COSMO BIO Co., LTD.

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル

URL:<http://www.cosmobio.co.jp/>

●営業部(お問合せ)

TEL:(03) 5632-9610

TEL:(03) 5632-9620

●札幌事業部(技術的なお問合せ)

TEL:(03) 5632-9744

FAX:(0134) 61-2295