

Human Type V Collagen

Cat. No. Y-4

Updated on Jan. 5th, 2024

【 I 】 Product Information

This is a purified human Type V Collagen.

Storage : -20°C

Description	Quantity	Concentration	Label	Form
Type V Collagen (Human)	150 µg	150 µg/mL	Unlabeled	5 mM Acetic acid

【 II 】 References

- [1] Niyibizi, C. et al., Human placenta type V collagens. Evidence for the existence of an α 1(V) α 2(V) α (3V) collagen molecule. J. Biol. Chem., 259, 14170-14174, 1984.

【 III 】 Related products

Description	Cat. No.	Quantity
Human Type I Collagen	Y-1	1 MG
Human Type III Collagen	Y-2	1 MG
Human Type IV Collagen	Y-3	150 µG
Human Type V Collagen	Y-4	150 µG
Human Type VI Collagen	Y-5	10 µG
Recombinant hTIMP-1	Y-8	100 µG
Recombinant hTIMP-2	Y-9	100 µG
Recombinant bTIMP-1	Y-10	100 µG
Recombinant hTIMP-1	Y-11	10 µG
Recombinant hTIMP-2	Y-12	10 µG
Recombinant bTIMP-1	Y-13	10 µG

**COSMO BIO Co., LTD.**

【JAPAN】
 TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
 KOTO-KU. TOKYO 135-0016, JAPAN
 Phone: +81-3-5632-9610
 FAX: +81-3-5632-9619
 URL: <https://www.cosmobio.co.jp/>

**COSMO BIO USA**

【Outside Japan】
 2792 Loker Ave West, Suite 101
 Carlsbad, CA 92010, USA
 email: info@cosmobiousa.com
 URL: www.cosmobiousa.com
 Phone/FAX: (+1) 760-431-4600