

CD11a Antibody : B-B15 FITC conjugated

HUMAN CD11a FITC CONJUGATED

CATALOGUE N° : 852.511.010

CLONE : B-B15

ISOTYPE : Mouse IgG1

IMMUNISATION : Human PBL

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognises the α L integrin of LFA-1 complex, a 180 kDa protein

QUANTITY : 100 tests / 1 ml

PURIFICATION : Ion exchange chromatography


FORMAT : Phosphate-buffered saline with 1 % BSA and 0.1 % sodium azide

WORKING DILUTION : Use 10 μ l to label 10^6 cells or 100 μ l of whole blood

STORAGE : Stable at +2-8°C for 12 months. For longer storage, freeze aliquots

APPLICATION : Suitable for flow cytometry

REFERENCES : Galon, J. et al., J Immunol.,1996; 157(3): 1184-92.
Gannon, G. A. et al., Am J Physiol.1998; 275(6 Pt 2): R1725-34.
Hahn, J. H. et al., J Immunol.,1997; 159(5): 2250-8.
Takasaki, J. et al., J Leukoc Biol.,1996; 60(2): 174-80.



A typical staining pattern with the B-B15 monoclonal antibody of lymphocytes

 Sodium azide is a poisonous and hazardous substance, handle and dispose of properly

Version 3 – 02.09

For research use only