

CD11a Antibody : B-B15 PE conjugated

HUMAN CD11a PE CONJUGATED

CATALOGUE N° : 852.512.010

CLONE : B-B15

ISOTYPE : Mouse IgG1

IMMUNISATION : Human PBL

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognises the α L integrin of LFA-1 complex, a 180 kDa protein

QUANTITY : 100 tests

PURIFICATION : Ion exchange chromatography


FORMAT : Phosphate-buffered saline with 1% BSA and 0.1% sodium azide

WORKING DILUTION : Reconstitute with 1 ml deionised water
Use 10 μ l to label 10^6 cells or 100 μ l of whole blood

STORAGE : Stable at +2-8°C for six months after reconstitution. DO NOT FREEZE

APPLICATION : Suitable for flow cytometry

REFERENCES : Canfield, S. M. and A. Y. Khakoo, J Immunol., 1999;163(6): 3430-40.
Ohashi, Y. et al., Am J Physiol., 1997; 272(6 Pt 2): H2807-14.



A typical staining pattern with the B-B15 monoclonal antibody of lymphocytes

 Sodium azide is a poisonous and hazardous substance, handle and dispose of properly

Version 3 – 02.09

For research use only