

CD34 Antibody : B-C34 PE conjugated

HUMAN CLASS I CD34 PE CONJUGATED

CATALOGUE N° : 854.592.010

CLONE : B-C34

ISOTYPE : Mouse IgG1

IMMUNISATION : KG-1a and TF1 cell lines

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognises the the Gp105-120 antigen, a 105-120 kDa protein

WORKSHOP : V

QUANTITY : 100 tests / 1 ml


PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline with 5 % BSA and 0.1 % sodium azide

WORKING DILUTION : Use 10 µl to label 10⁶ cells or 100 µl of whole blood

STORAGE : Stable at +2-8°C for 12 months
For longer storage, freeze aliquots

APPLICATION : Suitable for flow cytometry

 *Sodium azide is a poisonous and hazardous substance, handle and dispose of properly*

Version 2 - 11.06

For research use only