

CD63 Antibody : TS63 Azide free
HUMAN CD63 AZIDE FREE

CATALOGUE N° : 857.770.000

CLONE : TS63

ISOTYPE : IgG1

IMMUNISATION : Jurkat and HEL cell lines

HYBRIDOMA : Myeloma X63/AG8. x Balb/c spleen cells

SPECIFICITY : Recognises the LIMP antigen, a 30-60 kda (smear) protein

QUANTITY : 200 µg / 200 µl

PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

STORAGE : Stable at +2-8°C for 12 months
For longer storage, freeze aliquots

APPLICATION : Suitable for flow cytometry , immunoprecipitation and Western-Blot

REFERENCES : Charrin, S. et al., Eur J Immunol.,2003; 33(9): 2479-89.
Charrin S. et al., FEBS Lett., 2002; 516: 139-144
Charrin S. et al., J. Biol. Chem., 2001; 276: 14329-14337
Sordat I. et al., Lab. Invest., 2002; 82: 1715-1724

Version 2 - 01.08

For research use only