



Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Bovine Fibrinogen, Plasminogen Depleted

Concentration: 44.17 mg/ml

Aliquots: 1 x 1.00 g

Volume/aliquot: 22.640 ml

Activity: 100.00 % Clottable

Total protein: 1.00 g

Buffer: 20 mM sodium citrate-HCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 10/2006

PO: N/A

The Bovine Fibrinogen, Plasminogen depleted is >95% pure by SDS-PAGE and is greater than 95% clottable. The plasminogen was depleted using a lysine-sepharose column. Protein must be thawed in a 37 C water bath, aliquot and freeze the unused portion at <-60oC. All bovine material is from U.S. cows < 30 months of age.

Thank You,
Quality Assurance

Sample