



Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human Fibrinogen Plasminogen Depleted

Concentration: 52.72 mg/ml

Aliquots: 1 x 1.00 g

Volume/aliquot: 18.968 ml

Activity: 100.00 % Clottable

Total protein: 1.00 g

Buffer: 20 mM sodium citrate-HCl/ pH 7.4

Storage: < -60 °C

Expiration Date: 02/2007

PO: N/A

The Human Fibrinogen is homogeneous as judged on a 4-20% gradient gel and is greater than 95% clottable. The plasminogen was depleted using affinity chromatography. We suggest thawing this protein in a 37°C water bath, aliquoting and refreezing the unused portion. The above protein was purified from Human plasma that was tested and found negative for HIV-1, HIV-2, Hepatitis B Surface antigen and HCV. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample