



Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Human Antithrombin

Concentration: 1.89 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.529 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 20 mM Tris-HCl/ 0.1 M sodium
citrate/ 0.15 M NaCl/ pH 8.3

Storage: < -60 °C

Expiration Date: 05/2007

PO: N/A

The Human Antithrombin is a single band on SDS-PAGE and shows no reduction upon incubation with 2-mercaptoethanol. The above protein was purified from Human plasma that was tested and found negative for HIV-1, HIV-2, Hepatitis B Surface antigen and HCV. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample