



Androstenedione-3-CMO-HRP

Product type -
Form Liquid
Volume 1 ml (1000TEST)

Application notes ELISA
Recommended use together
Antibody : Anti Androstenedione-3-CMO-BSA (Cat.no. FKA-138E)
Recommended dilutions
ELISA : 50 µl/Well use

Storage Optimal dilutions/concentrations should be determined by the end user.
Maintain at 4°C-8°C in undiluted aliquots.
Avoid freeze/thaw cycles

References The Cross reactivities of anti Androstenedione-3-CMO (FKA-138E)
by ELISA with Androstenedione-3-CMO-HRP (FKA-137)

O.D. (492/600nm) of B0 about 0.7
Androstenedione dose of B0/2 about 5 ng/ml

Steroid

Cross Reaction

Androstenedion	100 %
5Alpha-Androstanedione	35.00 %
DHA	4.50 %
Androstenetrione	1.50 %
Progesteropne	1.20 %
Testosterone	6.00 %
Cortisol	0.32 %
Corticosterone	0.20 %
Estradiol	0.01 %



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619