



Pregnenolone-3-Succinate-HRP

Product type

Form

Liquid

Volume

1 ml (1000TEST)

Concentration

Application notes

ELISA

Recommended use

Recommended dilutions

ELISA : 50 µl/Well use

Storage

Optimal dilutions/concentrations should be determined by the end user.

Maintain at 4°C-8°C in undiluted aliquots. Avoid freeze/thaw cycles

References

*The Cross reactivities of anti Pregnenolone-3-Succinate (FKA-316E) by ELISA with 11-Oxo-Testosterone-3-CMO-HRP (FKA-315)

O.D. (492nm/600nm) of B₀ : 0.8

pregnenolone dose of B₀/2 : 1 ng/well

Steroid

Cross Reaction

Pregnenolone(P5)	100%
Pregnenolone-3- Succinate(P5-3-S)	37.00%
Progesterone(P4)	6.00%
17 alpha-Hydroxy-progesterone	1.00%
17 alpha-Hydroxy-progesterone	0.30%
DOC	0.08%
Cortisol(F)	0.01%
Cortisone(E)	0.02%
4-androstenedione	0.02%
Testosterone	0.02%



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619