



Aldosterone-3-CMO-HRP

Product type

Form

Liquid

Volume

1 ml (1000TEST)

Concentration

Application notes

ELISA

Recommended use

Recommended dilutions

ELISA : 50 µl/Well use

Storage

Optimal dilutions/concentrations should be determined by the end user.

Maintain at 4°C-8°C in undiluted aliquots. Avoid freeze/thaw cycles

References

*The Cross reactivities of anti Aldosterone-3-CMO-BSA (FKA-428E) by ELISA with Aldosterone-3-CMO-HRP (FKA-427).

O.D. (492nm/600nm) of B₀ : 0.6-1.0

Aldosterone dose of B₀/2 : 0.7-1.0 ng/ml (35-50 pg/well)

Steroid	Cross Reaction
Aldosterone	100%
Corticosterone (Compd B)	0.01%
DOC	0.02%
11-Deoxycortisol (compd S)	0.006%
Progesterone	0.00%
Cortisol (Compd F)	0.002%
17alpha-Hydroxyprogesterone	0.002%
Androstenedione	0.048%
Testosterone	0.002%
Dehydroepiandrosterone	0.010%
Estradiol	<0.001%



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp/index_e.asp e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619