



FD Section Storage Solution™

Product Name	Catalog #	Size
FD Section Storage Solution™	PC101	250 ml

The FD Section Storage Solution™ is a specifically formulated buffer solution that stays liquid at -20°C. It is ideal for the prolonged storage of free-floating fixed tissue sections. This solution has been proven to preserve the best morphology and antigenicity for cryostat and vibratome sections. Sections kept in this solution may be stored at -20°C for at least 5 years before processing for immunohistochemistry or routine histology.

References using FD Section Storage Solution™ (PC101):

1. Eid T, Lee TSW, Thomas MJ, Amiry-Moghaddam M, Bjørnsen LP, Spencer DD, Agre P, Ottersen OP, and de Lanerolle NC. (2005) Loss of perivascular aquaporin 4 may underlie deficient water and K⁺ homeostasis in the human epileptogenic hippocampus. *Proc. Natl. Acad. Sci. USA.* 102:1193-8.