

SPECIFICATION SHEET

Description:	Anti PSF1 monoclonal		
Catalog #:	GS001-100, GS001-50		
Immunogen:	mouse PSF1	Host:	Rat
Clone#:	aho57.2		
Subclass:	IgG		
Concentration:	1.0mg/ml		
Cross Reactivity:	PSF1, human 100%		
Form:	Purified by Protein G / Sep4B, supplied in phosphate buffered saline containing 0.09% NaN ₃		
Storage:	store at -20°C.		
Applications:	Suitable for use in immunohistochemistry, Immunocytochemistry and Western blotting.		

Additional Remarks:

Product without NaN₃ is also available on your request.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

11/06/01



COSMO BIO CO., LTD.

2-20, Toyo 2-Chome, Koto-ku, Tokyo 135-0016, JAPAN

T: +81-3-5632-9617

日本のお客様

F: +81-3-5632-9618

T: 03-5632-9610

E: export@cosmobio.co.jp

F: 03-5632-9619

W: www.cosmobio.com

W: www.cosmobio.co.jp