

SPECIFICATION SHEET

Description: Anti SLD5 monoclonal

Catalog #: GS005-100, GS005-50

Immunogen: mouse SLD5

Host: Rat

Clone#: AFS3-3-4

Subclass: IgG

Concentration: 1.0mg/ml

Cross Reactivity: SLD5, human 100%

Form: Purified by Protein G / Sep4B, supplied in phosphate buffered saline containing 0.09% NaN₃

Storage: store at -20°C.

Applications: Suitable for use in immunohistochemistry, Immunocytochemistry and Western blotting.

Additional Remarks:

Product without NaN₃ is also available on your request.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

11/06/01



COSMO BIO CO., LTD.

2-20, Toyo 2-Chome, Koto-ku, Tokyo 135-0016, JAPAN

T: +81-3-5632-9617

日本のお客様

F: +81-3-5632-9618

T: 03-5632-9610

E: export@cosmobio.co.jp

F: 03-5632-9619

W: www.cosmobio.com

W: www.cosmobio.co.jp