

## SPECIFICATION SHEET

**Description:** Anti SLD5 monoclonal

**Catalog #:** GS005-100, GS005-50

**Immunogen:** mouse SLD5

**Host:** Rat

**Clone#:** AFS3-3-4

**Subclass:** IgG

**Concentration:** 1.0mg/ml

**Cross Reactivity:** SLD5, human 100%

**Form:** Purified by Protein G / Sep4B, supplied in phosphate buffered saline containing 0.09% NaN<sub>3</sub>

**Storage:** store at -20°C.

**Applications:** Suitable for use in immunohistochemistry, Immunocytochemistry and Western blotting.

**Additional Remarks:**

Product without NaN<sub>3</sub> is also available on your request.

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**

11/06/01



**COSMO BIO CO., LTD.**

2-20, Toyo 2-Chome, Koto-ku, Tokyo 135-0016, JAPAN

T: +81-3-5632-9617

日本のお客様

F: +81-3-5632-9618

T: 03-5632-9610

E: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

F: 03-5632-9619

W: [www.cosmobio.com](http://www.cosmobio.com)

W: [www.cosmobio.co.jp](http://www.cosmobio.co.jp)