



仕 様 書

抗ヒト CEA モノクローナル抗体 (マウス)

コード 2AHCMC1

クローン番号 Hyb- 3519

由来 ヒト癌組織由来 Carchinoembryonic Antigen (CEA) を免疫した BALB/c マウスの脾臓細胞とマウスミエローマ細胞 (NS-1) を融合して得たハイブリドーマを BALB/c マウスに接種して得られた腹水

特異性 CEA+, NCA+, NCA II+の反応性を示す

精製 硫酸アンモニウム塩析精製

Ig Class IgG₁

溶媒 0.1% アジ化ナトリウムを含む生理的食塩水溶液

蛋白濃度 0.5 mg/mL

容量 0.2 mL

貯法及び有効期間 2~8°C 3ヶ月
-20°C 3年間 (凍結融解の繰り返しはさけてください)

用途 ELISA

参考文献

販売元



人と科学のステキな未来へ
コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
URL : <http://www.cosmobio.co.jp/>

● 営業部 (お問い合わせ)
TEL : (03) 5632-9610 FAX : (03) 5632-9619
TEL : (03) 5632-9620

製造元



〒112-0004 東京都文京区後楽1-1-10
Tel:03-3814-4081 Fax:03-3814-5957
e-mail:info@tokumen.co.jp
<http://www.tokumen.co.jp>

2AHCMC10806-02

Specification Sheet

Monoclonal Antibody Against Human CEA

Code 2AHCMC1

Clone No. Hyb- 3519

Manufacturing Process Hybridoma obtained by fusing mouse myeloma cells (NS-1) with BALB/c mouse spleen cells immunized with carcinoembryonic antigen (CEA) derived from human cancer tissue is proliferated in the BALB/c mouse abdominal cavity and ascites is collected.

Specificity The reactivity is CEA+, NCA+ and NCA II +.

Purification Ammonium sulfate precipitation.

Ig Class IgG₁

Form Saline (containing 0.1% NaN₃)

Protein Concentration 0.5 mg/mL

Package 0.2 mL

Storage 2-8°C (Below -20°C for long term storage)

Application ELISA

Reference

Distributor



COSMO BIO USA.

2792 Loker Avenue West, Suite 101, Carlsbad,
CA 92010 USA
Phone: +1 760-431-4600 / Fax: +1 760-431-4604
Email: info@cosmobiousa.com
URL: <https://www.cosmobiousa.com/>

Manufacturer



INSTITUTE OF IMMUNOLOGY CO., LTD.

1-1-10, Koraku, Bunkyo-ku, Tokyo, JAPAN 112-0004
Tel : +81-3-3814-4081 Fax : +81-3-3814-5957
e-mail : info@tokumen.co.jp
<http://www.tokumen.co.jp/english/>

2AHCMC10806-02

This product is for *in vitro* research use only and is not intended for use humans or animals.