



仕 様 書

THBR1-AS

コード	2ATHBR1
内容	マウス IgG ₁
由来	マウス腹水
性状	マウス腹水原液
防腐剤	0.1% アジ化ナトリウム
用途	マウス抗体を用いたイムノアッセイにおける、マウスに対する異好性抗体（HAMA:Human Anti-Mouse Antibody）による非特異反応を抑制します。 サンドイッチ測定系における標準添加濃度：0.1～1% 凝集法への応用も可能です。あらかじめ各測定系での最適条件を検討してください。
容量	1 mL
貯法及び有効期間	2～8℃ 3ヶ月 -20℃ 3年間 （凍結融解の繰り返しはさけてください）
備考	HBs抗原測定系には適しません。

販売元



人と科学のステキな未来へ
コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
URL: <http://www.cosmobio.co.jp/>

● 営業部（お問い合わせ）
TEL: (03) 5632-9610 FAX: (03) 5632-9619
TEL: (03) 5632-9620

製造元



株式会社 特殊免疫研究所

〒112-0004 東京都文京区後楽1-1-10
Tel: 03-3814-4081 Fax: 03-3814-5957
e-mail: info@tokumen.co.jp
<http://www.tokumen.co.jp>

2ATHBR11404-05

本製品は研究用にのみ使用できます。人やその他の動物の診断・治療に使用することはできません。

THBR1-AS (HAMA Blocker)

Code	2ATHBR1
Description	Mouse IgG ₁
Source	Mouse ascites
Purification	N/A
Preservative	0.1% NaN ₃
Application	This product inhibits nonspecific reaction caused by human anti-mouse antibody (HAMA) in immunoassays (ELISA, agglutination assay, and others) using mouse monoclonal antibodies. The recommended concentration is 0.1-1% for ELISA, however, the optimum condition should be determined for each assay.
Package	1 mL
Storage	2-8 °C (Below -20 °C for long term storage) Avoid multiple freeze and thaw cycles.
Note	This product is NOT appropriate for HBs antigen assays.

Distributor



2792 Loker Avenue West, Suite 101, Carlsbad,
CA 92010 USA
Phone: +1 760-431-4600 / Fax: +1 760-431-4604
Email: info@cosmobiousa.com
URL: <https://www.cosmobiousa.com/>

Manufacturer



INSTITUTE OF IMMUNOLOGY Co., LTD.

1-1-10, Koraku, Bunkyo-ku, Tokyo, JAPAN 112-0004
Tel : +81-3-3814-4081 Fax : +81-3-3814-5957
e-mail : info@tokumen.co.jp
<http://www.tokumen.co.jp/english/>

2ATHBR11404-05

This product is for *in vitro* research use only, and is not intended to use for diagnosis or treatment for humans nor animals.