

polyclonal antibody

ISL #126-3794

## Prion Protein

**Product code** 126-3794

**Region:**

**Alternative name**

**Additional info**

**Pack Size** 1 ml

**Sequence**

**Immunogen** KLH conjugate

**Host/Source** sheep

**Clone** n/a

**Subclass** n/a

**Format** purified Ig

**Preparation** prepared by Ammonium Sulphate fractionation

**Buffer System** PBS buffer, pH 7.2, 0.1% Sodium Azide

**State** Liquid

**Assay** ELISA 1/25 - 1/100

(BLOT?) to be determined

(IHC?) to be determined

**ELISA** This antibody is suitable for use in a simple ELISA, using 1 µg of unconjugated antigen coated per well, with alkaline phosphatase conjugated anti-sheep for detection. All incubations at 18°C using pNPP substrate, with OD measurements at 405 nm.

**End point** The end point titre for this antibody (in the above mentioned assay) is >1/1000.

**Storage** Aliquot and store at -20°C

*For in-vitro testing only. Not for diagnostic or in-vivo use*

*ISL (Immune Systems Ltd.) email@IS-L.com www.IS-L.com*

*Tel +44 1803 526 556 Fax +44 1803 526 776*