



COSMO BIO Co., LTD.
Inspiration for Life Science

K0602-M Anti DYKDDDDK Monoclonal Antibody (Clone No. 2H8)

Code No. K0602-M
Category Others
Target DYKDDDDK (N-terminal)
Type Monoclonal Antibody
Concentration 0.50mg/ml
Contents (Volume) 200 µg (400 µL/vial)
Gene ID -
Primary Source -
Synonyms -
Immunogen DYKDDDDK peptide expressing cell line


Raised in Mouse(BALB/c)
Myeloma SP2
Clone number 2H8
Purification ProteinG
Source Serum-free medium
Isotype IgG2a, κ
Cross Reactivity -
Label Unlabeled
Buffer PBS [containing 2% Block Ace as a stabilizer, 0.1% Proclin as a bacteriostat]
Storage Store below -20°C. Once thawed, store at 4°C. Repeated freeze-thaw cycles should be avoided.

Application WB, IP, ICC, IF, FCM (N-term only)

Recommended Antibody Dilutions


ELISA	WB	ICC	IF
Not Tested	0.01-1.0	Not Tested	0.1-1.0
IP	FCM	IF	Neutralization
0.1-1.0	0.1-10	0.25-5.0	—

(µg/mL)




ICC

Sample: CHO-K1/DYKDDDDK_hPGR Cell
Anti DYKDDDDK MAb (2H8) 0.1 µg/ml



FCM

Sample: DYKDDDDK_mBLT1 transfected CHO-K1
green: Standard Anti DYKDDDDK Mab 5 µg/ml
purple: Anti DYKDDDDK MAb (2H8) 5 µg/ml
gray: Control IgG



WB

Sample: CHO-K1/DYKDDDDK_hPGR Cell lysate 10ug
① Anti DYKDDDDK MAb (2H8) 0.1 µg/ml
② Standard Anti DYKDDDDK Mab 0.1 µg/ml
③ Control IgG 0.1 µg/ml

This antibody was prepared by Dr. Yokomizo T.
Department of Medical Biochemistry Department of Basic Medicine
Faculty of Medical Sciences Kyusyu university, Fukuoka, Japan

UniProt Summary

Reference

Manufactured by TransGenic Inc.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp/index_e.asp E-mail: export@cosmobio.co.jp

[Outside Japan]

[国内連絡先]

Phone : +81-3-5632-9617

Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619