



Certificate of Analysis

Product: Porcine Serum
Membrane Filtered

Lot #: E80605T

Catalog #: 100-115

Expiration Date: May 2004

Test/Method	Specification	Units	Results
pH		7.8	
Total Protein		7.95 mg%	
Hemoglobin		30.7 mg%	
Osmolality		347	
Sterility		Pass	

Product manufactured from raw material of U.S. Origin.

A stylized, handwritten signature in black ink, likely belonging to a representative of the Quality Control Department.

Quality Control Department



1301 East Beamer Street

Woodland, California 95776

(530) 668-3636

FAX (530) 668-3630

Porcine Serum

Sterile filtered

Mycoplasma tested and virus screened

Store frozen (-10°C to -20°C)

Porcine serum is collected from healthy adult animals.