

## Product Datasheet

[www.leinco.com](http://www.leinco.com)



### Mouse IgG2a Purified Low Endotoxin Functional Formulation Isotype Control

**Prod. No.:** I-118  
**Conc.:** 0.5 mg/ml  
**Pkg. Size:** Custom, 50 µg, 0.1 mg, 0.5 mg

#### Description

Specificity:

This IgG2a isotype control antibody has been tested against selected species' cells and tissues to assure minimal cross reactivity.

#### Host Species

Mouse

#### Formulation

This antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (PBS) pH 7.2, 150 mM NaCl with no carrier protein, potassium or preservatives added.

#### Purity

>98% by SDS Page and HPLC

#### Endotoxin

<0.5 EU/mg

#### Storage and Stability

This antibody is stable for at least one week when stored sterile at 2-8°C. For long term storage aseptically aliquot in working volumes without diluting and store at -20°C **in a manual defrost freezer. Avoid Repeated Freeze Thaw Cycles.**

#### Applications and Suggested Working Dilutions

This isotype control antibody should be used at the same concentration as the primary antibody.

#### Country of Origin

USA

#### Related Products



*Products are for research use only. Not for use in diagnostic or therapeutic procedures.*

359 Consort Drive | St. Louis, MO 63011 | (800) 538-1145 | (636) 527-5545 fax