



DNP-BSA

Form	Lyophilized protein powder 2,4-dinitrophenylated Bovine Serum Albumin
Volume	Protein: 20mg (As average of molecular weight: 69,000)
Storage	Store below 4°C (below -20°C for prolonged storage or after dissolved)
Concentration	Add distilled water. It will become 5mM phosphate buffered saline [PBS] , 0.9% NaCl solution.
Note	This product is lyophilized product containing totally 20mg proteins per vial. Therefore, conjugation ratio DNP groups/molecule, and volume of distilled water to reconstitute is different by production lots. For detailed information, visit the website below https://www.cosmobio.com/products/lsl_lg0017_1.html

Allergic antigen

LG-0009	DNP-Ascaris
LG-0533	DNP-Mite-Df
LG-0544	DNP-Mite-Dp
LG-0528	DNP-Cedar Pollen-Cj
LG-0017	DNP-BSA
LG-3017	DNP-BSA (30)
LG-0024	DNP-OA (Ovalbumin)
LG-0047	DNP-Casein
LG-0067	DNP-LG (BETA-Lactoglobulin)
LG-0089	DNP-KLH
LG-0317	BPO-BSA
LG-0577	BPO-BGG

LG-1109	TNP-Ascaris
LG-1117	TNP-BSA
LG-2334	Mite-Df Feces AG
LG-2444	Mite-Dp Feces AG
LG-5009	Ascaris Extract
LG-5339	Mite Extract-Df
LG-5449	Mite Extract-Dp
LG-5559	Mite Extract-Tp
LG-5229	Cedar Pollen Extract
LG-5280	Cedar Pollen Extract-Cj
LG-5779	Hinoki Cypress, Pollen Crude Extract
LG-5889	Ragweed Pollen Extract

Adjuvant

LG-6000	ALUM
---------	------

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

DNP-BSA

I. 内容

lot No. 518110

性 状	2,4-dinitrophenylated Bovine Serum Albumin 凍結乾燥粉末
規 格	蛋白質量：20mg DNP基数：1分子あたり8.8個（平均分子量69,000 としての計算値） 1. 2mlの精製水を加えることで5mM硼酸緩衝溶液 （pH8.0）0.9%NaCl溶液になります。


II. 保存上の注意

4℃以下（長期間又は溶液の場合は-20℃以下）で保存して下さい。

製造元

総発売元

 株式会社 エル・エス・エル

 コスモ・バイオ株式会社

新製品案内

LG-2334	Mite-Df Feces AG (ダニ-Df由来 排泄物抗原)
LG-2444	Mite-Dp Feces AG (ダニ-Dp由来 排泄物抗原)
LG-5999	Wormwood Pollen Extract
LG-7000	FTPP Suspension (百日咳死菌懸濁液)
LG-8000	FTMB Suspension (結核死菌懸濁液)
LB-7111	Anti (MAP) Mite Der fl
LB-7103	Anti (MAP) Mite Group 1
LB-7204	Anti (MAP) Mite Group 2