



DNP-LG (BETA-Lactoglobulin)

Form	Lyophilized protein powder 2,4-dinitrophenylated Bovine β -Lactoglobulins A&B
Volume	Protein: 30mg (As average of molecular weight: 18,400)
Storage	Store below 4°C (below -20°C for prolonged storage or after dissolved)
Concentration	Add distilled water. It will become 5mM phosphate buffered saline [PBS] , 0.9% NaCl solution.
Note	This product is lyophilized product containing totally 30mg proteins per vial. Therefore, conjugation ratio DNP groups/molecule, and volume of distilled water to reconstitute is different by production lots. For detailed information, visit the website below https://www.cosmobio.com/products/ls_lg0067_1.html

Allergic antigen

LG-0009	DNP-Ascaris
LG-0533	DNP-Mite-Df
LG-0544	DNP-Mite-Dp
LG-0528	DNP-Cedar Pollen-Cj
LG-0017	DNP-BSA
LG-3017	DNP-BSA (30)
LG-0024	DNP-OA (Ovalbumin)
LG-0047	DNP-Casein
LG-0067	DNP-LG (BETA-Lactoglobulin)
LG-0089	DNP-KLH
LG-0317	BPO-BSA
LG-0577	BPO-BGG

LG-1109	TNP-Ascaris
LG-1117	TNP-BSA
LG-2334	Mite-Df Feces AG
LG-2444	Mite-Dp Feces AG
LG-5009	Ascaris Extract
LG-5339	Mite Extract-Df
LG-5449	Mite Extract-Dp
LG-5559	Mite Extract-Tp
LG-5229	Cedar Pollen Extract
LG-5280	Cedar Pollen Extract-Cj
LG-5779	Hinoki Cypress, Pollen Crude Extract
LG-5889	Ragweed Pollen Extract

Adjuvant

LG-6000	ALUM
---------	------

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

LSL

カタログNo. LG-0067

DNP-LG (β -Lactoglobulin)

I. 内容

lot No. 892012

性 状	2,4-dinitrophenylated Bovine β -Lactoglobulins A&B 凍結乾燥粉末
規 格	蛋白質量: 30mg DNP基数: 1分子あたり 4.9個 (平均分子量 18,400 としての計算値) 1. 5mlの精製水を加えることで5mM硼酸緩衝溶液 (pH 8.0) 0.9% NaCl溶液になります。

II. 保存上の注意

4℃以下 (長期間又は溶液の場合は-20℃以下) で保存して下さい。

製造元

総発売元



株式会社 エル・エス・エル

コスモ・バイオ株式会社

関連製品

- ☆Ascaris Extract (アスカリス粗抽出物の凍結乾燥粉末)
カタログNo. LG-5009 (LSL)
- ☆DNP-Ascaris
カタログNo. LG-0009 (LSL)
- ☆DNP-BSA
カタログNo. LG-0017 (LSL)
- ☆DNP-OA (Ovalbumin) 新製品
カタログNo. LG-0024 (LSL)
- ☆DNP-Casein 新製品
カタログNo. LG-0047 (LSL)
- ☆DNP-LG (β -Lactoglobulin) 新製品
カタログNo. LG-0067 (LSL)
- ☆アラム (水酸化アルミニウムゲル) 新製品
カタログNo. LG-6000 (LSL)
- ☆抗DNP-Ascaris
カタログNo. LB-9009 (LSL)