

**CD222**  
**Clone MEM-238**  
**Cat. no. MNS-MON 2064**

---

**Specificity**

CD222 antigen (human). MEM-238 recognizes an epitope between domains 2 and 5. Immunogen: *Vaccinia-M6PR/iGFIIR*

**Immunoglobulin type**

Mouse IgG<sub>1</sub>

**Use**

The antibody can be used for Flow cytometry.

**Presentation**

100 µg in solution (PBS containing 15 mM sodium azide, pH 7.4).

**Literature**

- Leucocyte Typing VII. In press.

**Storage**

Store antibody at 2-8°C until expiry date.

FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.