

MDR Sampler Pack

Cat. no. MNS-MON 9200

Specificity

JSB-1 (ab-1) has potential value for diagnostic detection of MDR cells in human tumor samples.

LRP-56 (ab-2) has potential value for detection of LRP-associated non-Pgp MDR in human tumor samples.

MRPm6 (ab-3) and MRPr1 (ab-4) have potential value for detection of MRP-related non-Pgp MDR in human tumor samples. For more information look at individual data sheets.

Immunoglobulin type

JSB-1 (ab-1) murine IgG1

LRP-56 (ab-2) murine IgG2b

MRPm6 (ab-3) murine IgG1

MRPr1 (ab-4) rat IgG2a

Use

Immunocytochemistry: acetone-fixed cell preparations.

Immunohistochemistry: acetone-fixed frozen sections and formalin fixed paraffin embedded section.

Flow cytometry and Western blotting (Ab-1 can render less satisfactory results).

Presentation

This sampler pack contains 0.25 ml quantities of the regular products MON 9011, MON 9016, MON 9017 and MON 9018. All individual quantities are sufficient for a minimum of 50 tests.

Storage

Store at 4°C for short term (3 months) and at -20°C for extended storage.

FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.