

MCS

MOSS, INC.
P.O. BOX 189
PASADENA, MD 21122
1-800-932-6677

CERTIFICATE OF ANALYSIS

TMB
(3,3',5,5'-TETRAMETHYLBENZIDINE)
(in a proprietary buffer with enhancer)

PRODUCT- #TMBM-100,-500, -1000, 4 LITER, 10 LITER, 20 LITER

Expiration Date: JAN 1999 (18 months from date of manufacture at room temperature)
JULY 1999 (if refrigerated upon receipt)

LOT# 07198971

COLOR & CLARITY: clear	<input checked="" type="checkbox"/> _X_pass <input type="checkbox"/> __fail
pH 4.9 +/- 0.03	<input type="checkbox"/> -X-pass <input type="checkbox"/> ---fail
STABILITY @ 37°C :	<input type="checkbox"/> -X-pass <input type="checkbox"/> ---fail
STABILITY @ 18°C to 26°C	<input type="checkbox"/> -X-pass <input type="checkbox"/> ---fail
STABILITY @ 4°C	<input type="checkbox"/> -X-pass <input type="checkbox"/> ---fail
QC Raw Materials and log Lot#'s	<input type="checkbox"/> -X.yes <input type="checkbox"/> ---no
QC of lot during production (scans)	<input type="checkbox"/> -X.yes <input type="checkbox"/> ---no
Absorbance check of final product	<input type="checkbox"/> -X.yes <input type="checkbox"/> ---no
Performance data per "Immunoblot"	<input type="checkbox"/> -X-pass <input type="checkbox"/> ---fail

STORAGE CONDITIONS: 18°C to 26°C
FOR EXTENDED EXPIRATION DATE REFRIGERATE 2°C to 8°C

No special shipping conditions or precautions.

Ron Telford

Mmgr QA/QC