


TECHNICAL DATA SHEET**DMG-PEG2000, 1,2-dimyristoyl-rac-glycero-3-polyethylene glycol-2000**

Catalog Number: PG1-DG-2k

Product Structure:**Product Specifications:**

- Composition: DMG-PEG2000, 1,2-dimyristoyl-rac-glycero-3 (DMG)-PEG2000.
- Color: White.
- Form: Semi-solid/Solid.
- Lipid type: C14 fatty acid.
- PEG type: PEG2000, ~45 PEG units.
- Purity: >95%.

Storage Conditions:

DMG PEG2000 should be stored at -20 °C or lower. Desiccate. Other properties not known.

This product is for research use only and is not intended for use in humans or for diagnostic use.